# Pure Brite Electric Toothbrush

## User Instructions

*\*Read these instructions carefully and keep them for future reference.*

### Product Description

The Pure Brite Electric Toothbrush is an electric toothbrush that uses the vibration of a rechargeable motor to produce a bristle motion of both back and forth oscillation and rotation to clean your teeth. The vibrations easily remove germs and impurities around the teeth helping to maintain oral hygiene. The circular head design allows it to reach all around the tooth and clean 360° preventing any missed spots.

The acoustic energy produced by the motor transmits via the bristles to the teeth and gum's surface helping to loosen the bond between the dental plaque and your teeth, removing the bacteria in the gum pocket and tooth surface.

The long lasting battery allows for up to 45 days between charging.

### Caution

The Pure Brite Electric Toothbrush is suitable for people over 6 years of age. Children between 3 years or older must be supervised by an adult until they develop good brushing habits. Children under 3 years of age should not use this product.

If you had oral or gum surgery within 2 months consult your doctor or dentist prior to using this product.

Consult your dentist if you experience any discomfort or pain or are bleeding excessively after using this product.

### Cleaning and Care Instructions

The Pure Brite Electric Toothbrush has an IP level of IP65. Meaning Dust level 6: completely prevent dust from entering. Waterproof 5: Rinse with water without any harm.

1. Do not use water or acidic liquid above 104°F (40°C) to clean any part of the charging unit including the base, cable and plug.
1. DO NOT IMMERSE THIS PRODUCT IN WATER.
2. Remove the head and clean the head of the brush after every use with warm soapy water and rinse thoroughly.
  - a. If your toothpaste contains peroxide, sodium bicarbonate (baking soda) or other bicarbonates common in whitening toothpastes, be sure to remove all toothpaste from the head after every use to prevent plastic degradation.
3. Keep the charging interface clean and water-free when charging.
4. Replace brush head every 3 months for best performance. If the head is damaged replace right away.

## **Functions**

The Pure Brite Electric Toothbrush has 5 working modes:

1. Cleaning
2. Softening
3. Massage
4. Whitening
5. Nursing

## **How to use**

1. Press the 5S button to turn on, switch between modes and turn off the toothbrush.
2. The toothbrush will work for 2 minutes, pausing every 30 seconds to remind you to switch position. After 2 minutes of use, the toothbrush will turn off to prevent damage to your teeth and gums.

## **Warranty**

The Pure Brite Electric Toothbrush has a 3 year manufacturer defect warranty.

This product is only designed for cleaning teeth, gums and tongue. Use for any other purpose will void this warranty.

For use only with Pure Brite Brush heads, use with any other brand will void this warranty.

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Home Aide Diagnostics, Inc.