

Pr♥Health Mini

Talking Blood Pressure Monitor




Oscillometry results:

- mmHg or kPA
- Hypertension
- Hypotension
- Pulse (40~195) BPM
- 99 Memory Sets



Audible Results



2 Users



ID: 50632-0007-26
30bx/case

- Weight: ~0.5lb (excludes batteries)
- Size: 4.5x3.8x2.3 in (LxWxH)
- Cuff : 8.66-12.60 in (Circumference)



Home Aide Diagnostics, Inc.
1072 S Powerline Rd.
Deerfield Beach, FL 33442
1-800-915-0116 - www.homeaide.us
Made in China