



# PureAir

## Mini Nebulizer



- ◆ Portable and lightweight
- ◆ Extra fine mist delivery system
- ◆ 7ml medication cup
- ◆ Adult and Child Mask included
- ◆ Mouthpiece included
- ◆ AC power adapter included
- ◆ NDC# 50632-0007-55

Model	MCN-S600A
Product	Pure Air Mini Compressor
Power supply	Power adapter(AC100~240V 50H/60 DC 12V 1A)
Power consumption	≤ 11W
Compressor pressure Range	Approx.16 ~22 Psi/110 ~ 151 Kpa
Operation pressure Range	Approx.6 ~ 10 Psi/41 ~69Kpa
Liter flow Rang	Approx.2.5 ~ 5 lpm
Nebulization rate	≧ 0.20 ml/min
Particle size	Approx.(72.68±0.39)% within 0.5 to 5um range
MMAD	2.93±0.03um
Capacity of medication cup	Max. 7ml
Noise level	≤49dBA(acc.to DIN EN 13544-1 section 26)
Expected service life	400h
Operating temperature and humidity	10 ~40°C, 30 ~ 85 % RH, 800~1060hPa
Storage and delivery condition	-20~ 70°C, 20 ~ 75 % RH, 800~1060hPa
Safety level	Type B ClassII
Dimension(LxWxH)	Approx. 117x90x58mm
Weight	Approx.210g



**Home Aide Diagnostics, Inc.**

1072 S. Powerline Road \* Deerfield Beach, FL 33442

www.homeaide.us \* 1 (800) 915 - 0116