



PureAir

Mini Nebulizer



- ◆ Portable and lightweight
- ◆ Extra fine mist delivery system
- ◆ 7ml medication cup
- ◆ Adult and Child Mask included
- ◆ Mouthpiece included
- ◆ AC power adapter included
- ◆ NDC# 50632-0007-55

Product	DC Compressor Nebulizer
Power Supply	Power adapter (input AC 100~240V 50/60Hz 350mA, output 12V 1A)
Power Consumption	≤ 14W
Compressor Pressure Range	Approx.16~22 Psi/110~151 Kpa
Operation Pressure Range	Approx.6~10 Psi/41~69Kpa
Liter Flow Rang	Approx.2.5~5 lpm
Nebulization Rate	≧0.25 ml/min
Particle Size	Approx.(72.68±0.39)% within 0.5 to 5um range
MMAD	2.93±0.03um
Capacity of Medication Cup	Max. 7ml
Noise level	≤ 49.2dBA(acc.to DIN EN 13544-1 section 26)
Expected Service Life	400h
Dimension	Approx. L160X W130 X H110 mm
Weight	Approx.390g
Operating Temperature and Humidity	10~40°C, 30~85% RH, 800~1060hPa
Storage and Delivery Condition	-20~70°C, 20~75% RH, 800~1060hPa
Safety Level	Type B ClassII



Home Aide Diagnostics, Inc.

1072 S. Powerline Road * Deerfield Beach, FL 33442

www.homeaide.us * 1 (800) 915 - 0116